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Product version	DMA 35	DMA 35 Tag&Log	DMA 35 Ex	DMA 35 Ex Petrol
Measuring range	Density: 0 to 3 g/cm ³ Temperature: 0 °C to 40 °C (32 to 104 °F) Viscosity: 0 to 1000 mPa·s			
Accuracy	Density*: 0.001 g/cm ³ Temperature: 0.2 °C (0.4 °F)			
Repeatability	Density: 0.0005 g/cm ³ Temperature: 0.1 °C (0.2 °F)			
Resolution	Density: 0.0001 g/cm ³ Temperature: 0.1 °C (0.1 °F)			
Ambient temperature**	-10 °C to +50 °C (14 to 122 °F)			

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Meas. units	Options	Description
Density	Density	Density at the displayed measuring temperature in g/cm ³ or kg/m ³ .
	Density @ xx °C: α (g/cm ³ /K)	Density at the chosen reference temperature in g/cm ³ or kg/m ³ . The temperature influence is compensated by the set temperature coefficient α .
	Specific Gravity SG: SG Temp. (°C or °F) α (g/cm ³ /K)	Specific gravity is the density of the sample (at the chosen reference temperature) divided by the density of water (at the chosen reference temperature). The temperature influence is compensated by the set temperature coefficient α (g/cm ³ /K).
Alcohol	Alcohol % v/v @ 20 °C	Concentration of an ethanol/water mixture in % by volume at 20 °C. The measuring range is 0 to 100 % v/v.
	Alcohol % w/w	Concentration of a mixture of alcohol and water in % by weight. The measuring range is 0 to 100 % w/w.
	Alcohol US @ 60 °F (°Proof)	Degrees Proof at 60 °F. The measuring range is 0 to 200 °Proof.
API ^a	API Gravity A API Gravity B API Gravity D	API number for the product group referred to the reference temperature of 15 °C or 60 °F. Product group A: Crude oil Product group B: Fuels Product group D: Lubricants
	API SG A API SG B API SG D	Specific gravity for the product group referred to the reference temperature of 15 °C or 60 °F.
	API Density A API Density B API Density D	Density of the product group in g/cm ³ with regard to the reference temperature of 15 °C or 60 °F.

* Depending on the set temperature unit (°C or °F), the API values are automatically displayed @ 15 °C or @ 60 °F.

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V 1.79

Anton Paar

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Meas. units	Options	Description
Baumé	α (g/cm ³ /K)	Degrees Baumé at 60 °F. The measuring range is 0 to 100 °Baumé. Temperature influence is compensated by the set temperature coefficient α . For determining degree Baumé there are two different calculation methods depending on whether the density is above or below the density of water. Depending on the density of the measured liquid, DMA 35 automatically switches between the two calculation methods.
H2SO4	H2SO4 % w/w	Concentration of sulfuric acid or battery acid in % by weight. The measuring range is 0 to 70 % w/w.
	H2SO4 @ 20 °C	Density of sulfuric acid or battery acid at 20 °C. The measuring range is 0.8 to 2 g/cm ³ .
Sugar	Brix	Degrees Brix (sucrose concentration in % by weight). The measuring range is -10 to 85 °Brix.
	Extract (°Plato)	Degrees Plato. The measuring range is -10 to 85 °Plato.
Custom Functions	--	Optional custom functions. Contact your Anton Paar representative to obtain a custom function.
Period	--	Period value of the oscillator at the measuring temperature.
Raw data	--	Period value and resistance of the temperature sensor (only for service purposes).